FICHE TECHNIQUE

COSMÉTIQUE ET PARAMÉDICAL





Ref. A12592 WALL-MOUNTED ELBOW DISPENSER KING

Ideal for preventing the spread of viruses and bacteria thanks to a stainless steel elbow-operated lever.

√ Compatible with 1 L ref. A19253 hand sanitizer refill



SPECIFICATIONS

Dimensions: 14 cm Material: stainless steel

INSTRUCTIONS FOR USE

- 1. Mount the pump to the wall
- 2. Place the stopper on the bottle
- 3. Insert the bottle into the pump
- 4. Use your elbow to activate the lever, which presses on the stopper, and collect the product in your hand.

PACKAGING

Sales unit: 1 unit

Subpackaging: box of 20 units

Packaging: box of 100 units (5 boxes of 20 units)

REGULATIONS

FOR INDIVIDUAL AND PROFESSIONAL USE.

Our data sheets are prepared with current knowledge and are subject to modifications. Photos non-binding. 20/06/2023 Index 1.0





