



ALCOHOL AT 70°

A12831 - 50ML

Disinfectant spray for external use.

- Body disinfectant (on healthy skin)
- Surface disinfectant (small equipment)

STANDARDS	STRAINS	CONDITIONS
Bact. EN 1276	<i>Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Enterococcus hirae</i>	In 5 minutes
Bact. EN 13727		In 1 minute

● Application:

For antisepsis of healthy skin. Can also be used on small items such as tweezers, small scissors, etc.

● Directions for use:

Apply the disinfectant directly to healthy skin. Leave to act for 5 minutes.

PRECAUTIONS FOR USE:

Avoid contact with eyes.

● Areas of use:

- Hospitals
- Retirement homes, EHPAD
- Maternity homes, nurseries and schools
- Campsites
- Accommodation, hotels

● Packaging:

Packaging unit (U.C): 50mL aerosol

Sales unit: Box of 24 units

Palettization: 4608 units / 192 boxes

● Legislation:

GENERAL PUBLIC AND PROFESSIONAL USE

Our technical sheets have been produced according to the current state of our knowledge; some elements are subject to change. Non-contractual photo. Our technical data sheets have been produced based on our current knowledge; certain elements are subject to change. Non-contractual photo. 01/06/2026 Index 1